Electronic Instructions for Use

eIFU Indicator

STATUS OF THE INFORMATION

27.08.2024

TRADE NAME

neotivCare

DATE OF MANUFACTURE

2024

VERSION NUMBER

Version 1.7.1

MANUFACTURER

neotiv GmbH

Hegelstr. 19

39104 Magdeburg

www.neotiv.com

INSTRUCTIONS FOR USE

These instructions for use are intended to provide additional information and support the use of the neotivCare software for smartphones and tablets (app). Patients can download and install the application themselves on their own device. The use of neotivCare is safe and self-explanatory even without instructions for use. The instructions in this document can also be taken into account during use. If you have any questions, please contact the neotiv support team at support@neotiv-care.com.

I. GENERAL INFORMATION

INTENDED USE

The interactive software application "neotivCare" is intended to enable an assessment of cognitive performance to detect and monitor cognitive impairment in adults, and in particular to support the diagnosis of a mild cognitive impairment or an early Alzheimer's dementia in 60- to 80-year-olds.

Intended users perform the tests independently in their domestic environment. The report obtained is intended to be discussed with the supervising physician and support diagnostic and therapeutic decision.

HOW IT WORKS

neotivCare is designed for use by patients and is based on digital cognitive tests especially targeting memory functions. These tests are specifically designed for brain areas that are affected early by cognitive impairment and are related to an early stage of Alzheimer's disease. The mobile software application neotivCare (app) is therefore primarily intended to be used by patients aged 60-80 years to objectively assess cognitive impairment and deterioration especially in the area of memory functions. Adult patients outside the specified 60-80 age group can also use the app to assess their cognitive performance. The assessment of the cognitive performance of these age groups is only possible indirectly in relation to the reference group of 60–80-year-olds and must be interpreted with appropriate caution by the physician. The assessment of cognitive impairment is based on 3 different interactive, visual cognitive tests that are completed unsupervised, independent of the environment, ideally in the home environment, on a mobile device (smartphone or tablet). The patient can set the preferred time and a reminder for the upcoming tests and is informed by push notifications when new tests are available. Before the first test can be performed, training must be completed within the app. During the tests, scenes and objects or differences in scenarios must be identified. Answer options are displayed and entered by tapping on the screen of the smartphone or tablet. The tapped area is registered in the app by measuring the incorrect and correct entries (touches) of the patient and thus enables an objective recording of the measurement results. A composite total score is calculated from the results of the test sessions, based on at least 2 out of 4 valid results of the test sessions. The composite score is compared with the scores of the reference group or a previous assessment of the same person. In addition, after completing each test, the patient is asked to answer a self-assessment about the test situation, the results obtained and the external circumstances. In the app, the patient can create a printable letter of findings that includes the score evaluation and self-report answers.

PERFORMANCE PARAMETERS

neotivCare tests mechanisms of cognition (especially memory) that are neuroanatomically validated to detect cognitive impairment related to specific brain dysfunction.

The comparison of the test results with the reference group allows to distinguish between unremarkable cognitive performance and cognitive impairment typical of mild cognitive impairment (amnestic variant) associated with Alzheimer's disease in patients aged 60-80 years with a sensitivity of 79% and a specificity of 75%. The sensitivity and specificity are based on the identification of cognitive impairment e.g. by the preclinical Alzheimer's cognitive composite score (PACC5).

The cognitive performance of patients aged 18-59 years can be monitored throughout the prescription period. Changes in cognitive impairment in these patients can only be viewed in comparison to their own previous test results.

PREREQUISITE FOR USE AND TECHNICAL REQUIREMENTS FOR DEVICES TO BE COMBINED

neotivCare should only be used by patients, who:

1. are confident in handling of smartphones/tablets
2. have a stable wireless internet connection (wifi).
3. do not suffer from physical limitations that prevent the safe use of smartphones/tablets.
4. have sufficient language skills in German or English.

Access to the app is gained after installation on the end device by creating a user account with a user name and password as well as an email address to be verified. The tests can be started after entering a 16-digit activation code. This activation code is made available to the user, for example after submitting a prescription to his or her health insurance company.

The technical requirements for using the app are a mobile device, i.e. a smartphone or tablet with a touch screen and a compatible operating system (from Android version 9.0 (2018) or from iOS version 12.0 (2017)), as well as access to the internet via WLAN (WiFi) or mobile data networks.

TYPE AND DURATION OF USE

The planned duration of use of neotivCare is 12 weeks. The period of use can be extended by postponing the tests. Use can be terminated at any time. The initial use can be followed by two further 12-week periods of use.

FREQUENCY OF USE

During the period of use, weekly tests are carried out, which consist of a learning phase and a test phase that are separated in time.

MINIMUM AGE

18 years of age

SIDE EFFECTS

There are no known side effects, if you suspect any side effects, please report them to support@neotiv-care.com.

INTERACTIONS

There are no known interactions. If you notice any interactions, please report them at support@neotiv-care.com.

PAPER-BASED INSTRUCTIONS FOR USE

You will receive a complimentary paper version of the instructions for use within seven working days on request by e-mail at support@neotiv-care.com.

II. INFORMATION FOR PATIENTS

GENERAL INFORMATION

If you perceive a rapid deterioration in your cognitive abilities, we recommend that you consult a doctor immediately and, if necessary, discontinue the test with neotivCare.

Self-perceived memory problems and intense preoccupation with the subject of memory problems and dementia may in individual cases lead to increased anxiety, fear and depression. In this case, please contact your attending physician immediately.

SPECIAL INSTRUCTIONS FOR USE

Ensure that the environment in which the neotivCare app tests are performed is suitable for safe performance and use of the mobile device. Avoid using it in dangerous situations, for example while driving a motor vehicle.

A quiet and undisturbed environment also favours the attentive performance of the tests and leads to useful test results. Set a test time in the app that can be easily integrated into your everyday life, for example free time in the afternoon.

The use of mobile data networks may result in charges for which neotiv GmbH accepts no liability.

REPORTING OF SERIOUS INCIDENTS

If serious incidents occur in connection with the use of the medical device neotivCare, please contact our support team immediately at support@neotiv-care.com and report the incident to the Federal Ministry for Drugs and Medical Devices (BfArM). Serious incidents are defined as damage to health in connection with the use of a medical device that require medical treatment.

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Kurt-Georg-Kiesinger-Allee 3 53175 Bonn

Telefax: 0228/207 5300

You can also submit the notification online at:

<https://www2.bfarm.de/medprod/mpsv/>

DELETION OF USER ACCOUNT

If your period of use ends and no further activation code is used, your user account will be deleted after 3 years and all associated data anonymized.

You can delete your user account at any time in the app under the menu item "User account", "Personal settings", "Delete user account". This will immediately and irrevocably delete all data generated by you from your device and the server.

In both cases, you will no longer be able to log in with your access data to use the app after deleting your user account.

III. INFORMATION FOR PHYSICIANS

INDICATIONS

neotivCare is intended for the assessment of self-reported or suspected memory impairment in adults. The software is to be used in cases of suspected:

F00* Dementia in Alzheimer's disease (G30.-†).

F00.0* Dementia in Alzheimer's disease, with early onset (G30.0†)

F00.1* Dementia in Alzheimer's disease, with late onset (G30.1†)

F00.2* Dementia in Alzheimer's disease, atypical or mixed form (G30.8†)

F00.9* Dementia in Alzheimer's disease, unspecified (G30.9†)

F01 Vascular dementia

F01.0 Vascular dementia with acute onset

F01.1 Multi-infarct dementia

F01.2 Subcortical vascular dementia

F01.3 Mixed cortical and subcortical vascular dementia

F01.8 Other vascular dementia

F01.9 Vascular dementia, unspecified

F02.0* Frontotemporal dementia (G31.0)

F02.8* Lewy body dementia (G31.8) F06.7

Mild cognitive impairment

The ICD-10 classifies diagnoses primarily according to their aetiology. However, for some diagnoses it is necessary to additionally code manifestations that are not included in the code for the aetiology. This combined ordering principle is made possible by the so-called "cross-star system" of the ICD-10. First the aetiology (i.e. the primary code with a cross †) is coded, then the manifestation (= secondary code with an asterisk *).

CONTRAINDICATIONS


There are no known contraindications.

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